

## Product datasheet for **GA203813**

### Mouse **Cx3cl1** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	Cx3cl1 CRISPRa kit - CRISPR gene activation of mouse chemokine (C-X3-C motif) ligand 1
Symbol:	Cx3cl1
UniProt ID:	<u><a href="#">Q35188</a></u>
Synonyms:	AB030188; ABCD-3; AI848747; CX3C; Cxc3; D8Bwg0439e; Scyd1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)