

Product datasheet for **GA203837**

Mouse Sema4c activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Sema4c CRISPRa kit - CRISPR gene activation of mouse sema domain, immunoglobulin domain (Ig), transmembrane domain (TM) and short cytoplasmic domain, (semaphorin) 4C
Symbol:	Sema4c
UniProt ID:	Q64151
Synonyms:	AI426163; M-Sema F; Semacl1; Semaf; Semai; sema I
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)