

## Product datasheet for **GA203880**

### Mouse Shc3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Shc3 CRISPRa kit - CRISPR gene activation of mouse src homology 2 domain-containing transforming protein C3
Symbol:	Shc3
Components:	GA203880G1, Shc3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA203880G2, Shc3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA203880G3, Shc3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u><a href="#">Q61120</a></u>
Synonyms:	N-Shc; Rai; ShcC
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)