


Product datasheet for **GA203889**

Mouse Siah1a activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Siah1a CRISPRa kit - CRISPR gene activation of mouse siah E3 ubiquitin protein ligase 1A
Symbol:	Siah1a
Components:	GA203889G1, Siah1a gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA203889G2, Siah1a gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA203889G3, Siah1a gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P61092</u>
Synonyms:	AA982064; AI853500; Sinh1a
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)