

Product datasheet for **GA203890**

Mouse Siah1b activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Siah1b CRISPRa kit - CRISPR gene activation of mouse siah E3 ubiquitin protein ligase 1B
Symbol:	Siah1b
Components:	GA203890G1, Siah1b gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA203890G2, Siah1b gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA203890G3, Siah1b gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q06985
Synonyms:	AA960570; siah-1b; Sinh1b
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)