

Product datasheet for **GA203909**

Mouse Sin3a activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Sin3a CRISPRa kit - CRISPR gene activation of mouse transcriptional regulator, SIN3A (yeast)
Symbol:	Sin3a
UniProt ID:	Q60520
Synonyms:	AW553200; mKIAA4126; mSin3A; Sin3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)