

Product datasheet for **GA203944**

Mouse Slc20a1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Slc20a1 CRISPRa kit - CRISPR gene activation of mouse solute carrier family 20, member 1
Symbol:	Slc20a1
Components:	GA203944G1, Slc20a1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA203944G2, Slc20a1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA203944G3, Slc20a1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q61609
Synonyms:	AI607883; Glvr-1; Glvr1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)