


## Product datasheet for GA203945

### Mouse Slc20a2 activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | Slc20a2 CRISPRa kit - CRISPR gene activation of mouse solute carrier family 20, member 2 |
| Symbol:       | Slc20a2  |
| UniProt ID:   | <u><a href="https://uniprot.org/Q80UP8">Q80UP8</a></u>                                   |
| Synonyms:     | MolPit2; Pit-2; Pit2; Ram-1; Ram1  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                  |

#### Product images:



[View online >](#)