

Product datasheet for **GA203960**

Mouse Slc4a1ap activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Slc4a1ap CRISPRa kit - CRISPR gene activation of mouse solute carrier family 4 (anion exchanger), member 1, adaptor protein |
| Symbol: | Slc4a1ap |
| Components: | GA203960G1, Slc4a1ap gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA203960G2, Slc4a1ap gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA203960G3, Slc4a1ap gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077 |
| Synonyms: | kanadaptin |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online »](#)