

Product datasheet for **GA203962**

Mouse Slc4a3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Slc4a3 CRISPRa kit - CRISPR gene activation of mouse solute carrier family 4 (anion exchanger), member 3
Symbol:	Slc4a3
Components:	GA203962G1, Slc4a3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA203962G2, Slc4a3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA203962G3, Slc4a3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	P16283
Synonyms:	A930038D23Rik; Ae3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)