

Product datasheet for **GA203993**

Mouse *Sstr3* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Sstr3 CRISPRa kit - CRISPR gene activation of mouse somatostatin receptor 3
Symbol:	Sstr3
UniProt ID:	P30935
Synonyms:	Smstr-3; Smstr3; sst3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)