

Product datasheet for **GA203997**

Mouse Siglec1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Siglec1 CRISPRa kit - CRISPR gene activation of mouse sialic acid binding Ig-like lectin 1, sialoadhesin
Symbol:	Siglec1
UniProt ID:	Q62230
Synonyms:	Cd169; Siglec-1; Sn
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)