

Product datasheet for **GA204004**

Mouse Snap23 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Snap23 CRISPRa kit - CRISPR gene activation of mouse synaptosomal-associated protein 23 |
| Symbol: | Snap23 |
| UniProt ID: | O09044 |
| Synonyms: | 23kDa; AA408749; SNAP-23; Sndt; Syndet |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)