

Product datasheet for **GA204007**

Mouse Snrk activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Snrk CRISPRa kit - CRISPR gene activation of mouse SNF related kinase
Symbol:	Snrk
Components:	GA204007G1, Snrk gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA204007G2, Snrk gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA204007G3, Snrk gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8VDU5</u>
Synonyms:	2010012F07Rik; AI448042; AW547029; E030034B15; mKIAA0096; R74830
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)