

Product datasheet for **GA204020**

Mouse Son activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Son CRISPRa kit - CRISPR gene activation of mouse Son DNA binding protein
Symbol:	Son
Components:	GA204020G1, Son gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA204020G2, Son gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA204020G3, Son gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9QX47</u>
Synonyms:	2900011L12Rik; AA409051; AU067731; C81487; mKIAA1019; nrebp
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)