

Product datasheet for **GA204043**

Mouse Spa17 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Spa17 CRISPRa kit - CRISPR gene activation of mouse sperm autoantigenic protein 17
Symbol:	Spa17
Components:	GA204043G1, Spa17 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA204043G2, Spa17 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA204043G3, Spa17 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q62252</u>
Synonyms:	Sp17
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)