

Product datasheet for **GA204045**

Mouse Sp4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Sp4 CRISPRa kit - CRISPR gene activation of mouse trans-acting transcription factor 4
Symbol:	Sp4
Components:	GA204045G1, Sp4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA204045G2, Sp4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA204045G3, Sp4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	5730497N03Rik; HF-1b; HF1-b
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)