

Product datasheet for **GA204088**

Mouse **Sprr2a1** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Sprr2a1 CRISPRa kit - CRISPR gene activation of mouse small proline-rich protein 2A1
Symbol:	Sprr2a1
Components:	GA204088G1, Sprr2a1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA204088G2, Sprr2a1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA204088G3, Sprr2a1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q4KL71</u>
Synonyms:	Sprr2a; Sprr2a3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)