

Product datasheet for **GA204142**

Mouse Plk4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Plk4 CRISPRa kit - CRISPR gene activation of mouse polo like kinase 4
Symbol:	Plk4
UniProt ID:	Q64702
Synonyms:	1700028H20; AI385771; Sak; Stk18
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)