

Product datasheet for GA204155

Mouse Stx4a activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Stx4a CRISPRa kit – CRISPR gene activation of mouse syntaxin 4A (placental)
Symbol:	Stx4a
Components:	GA204155G1, Stx4a gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA204155G2, Stx4a gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA204155G3, Stx4a gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P70452</u>
Synonyms:	Stx4; Syn-4; Syn4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

