

## Product datasheet for **GA204175**

### Mouse Swap70 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Swap70 CRISPRa kit - CRISPR gene activation of mouse SWA-70 protein
Symbol:	Swap70
UniProt ID:	<a href="#">Q6A028</a>
Synonyms:	70kDa; AV235546
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)