


Product datasheet for **GA204191**

Mouse T activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	T CRISPRa kit - CRISPR gene activation of mouse brachyury, T-box transcription factor T
Symbol:	T
Components:	GA204191G1, T gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA204191G2, T gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA204191G3, T gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P20293</u>
Synonyms:	Bra; cou; D17Mit170; Low; Lr; me75; T1; Tl2; Tl3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)