

Product datasheet for **GA204193**

Mouse *Tac2* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tac2 CRISPRa kit - CRISPR gene activation of mouse tachykinin 2
Symbol:	Tac2
Components:	GA204193G1, Tac2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA204193G2, Tac2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA204193G3, Tac2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P55099</u>
Synonyms:	PPT-B; Tac3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

