

Product datasheet for **GA204195**

Mouse *Tacr1* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tacr1 CRISPRa kit - CRISPR gene activation of mouse tachykinin receptor 1
Symbol:	Tacr1
UniProt ID:	P30548
Synonyms:	Nk1r; Spr; Tac1r
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)