

Product datasheet for **GA204197**

Mouse **Tacr3** activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tacr3 CRISPRa kit - CRISPR gene activation of mouse tachykinin receptor 3
Symbol:	Tacr3
UniProt ID:	<u>P47937</u>
Synonyms:	Nk3r; Nkr; Tac3r
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)