

## Product datasheet for **GA204198**

### Mouse **Taf1a** activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Taf1a CRISPRa kit - CRISPR gene activation of mouse TATA-box binding protein associated factor, RNA polymerase I, A
Symbol:	Taf1a
UniProt ID:	<a href="#">P97357</a>
Synonyms:	mTAFI48; TAFI48
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)