

Product datasheet for **GA204199**

Mouse *Taf1b* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Taf1b CRISPRa kit - CRISPR gene activation of mouse TATA-box binding protein associated factor, RNA polymerase I, B
Symbol:	Taf1b
UniProt ID:	P97358
Synonyms:	4930408G01Rik; A230108M10Rik; p63; TAFI68; Tafi86
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)