

Product datasheet for **GA204201**

Mouse *Taf6* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Taf6 CRISPRa kit - CRISPR gene activation of mouse TATA-box binding protein associated factor 6
Symbol:	Taf6
UniProt ID:	Q62311
Synonyms:	80kDa; AW549759; p80; TAF(II)80; Taf2e; TAFII70
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)