

Product datasheet for **GA204207**

Mouse Tank activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Tank CRISPRa kit - CRISPR gene activation of mouse TRAF family member-associated Nf-kappa B activator |
| Symbol: | Tank |
| UniProt ID: | P70347 |
| Synonyms: | C86182; E430026L09Rik; I-TRAF |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:

