

Product datasheet for **GA204259**

Mouse **Tcp10b** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Tcp10b CRISPRa kit – CRISPR gene activation of mouse t-complex protein 10b
Symbol:	Tcp10b
Components:	GA204259G1, Tcp10b gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA204259G2, Tcp10b gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA204259G3, Tcp10b gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	D17Leh66B; D17Leh66ba; T66B-α; Tcp-10b; Tcp-10bt
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

