

Product datasheet for **GA204260**

Mouse **Tcp1l** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Tcp1l CRISPRa kit - CRISPR gene activation of mouse t-complex protein 1l
Symbol:	Tcp1l
Components:	GA204260G1, Tcp1l gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA204260G2, Tcp1l gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA204260G3, Tcp1l gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	D17Ken1; Tcp-1l
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

