

Product datasheet for **GA204282**

Mouse Tead3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tead3 CRISPRa kit - CRISPR gene activation of mouse TEA domain family member 3
Symbol:	Tead3
UniProt ID:	P70210
Synonyms:	DTEF-1; ETFR-1; Tcf13r2; TEAD-3; TEF-5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)