

Product datasheet for **GA204283**

Mouse Tead4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tead4 CRISPRa kit - CRISPR gene activation of mouse TEA domain family member 4
Symbol:	Tead4
UniProt ID:	Q62296
Synonyms:	ETFR-2; Etf2; FR-19; Rtef1; Tcf13r1; TEAD-4; TEF-3; Tef3; Tefr; Tefr1; Tefr1a
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)