

Product datasheet for **GA204284**

Mouse Alyref activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Alyref CRISPRa kit - CRISPR gene activation of mouse Aly/REF export factor
Symbol:	Alyref
UniProt ID:	O08583
Synonyms:	ALY; REF1; Refbp1; Thoc4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)