

Product datasheet for **GA204289**

Mouse Tek activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | Tek CRISPRa kit - CRISPR gene activation of mouse TEK receptor tyrosine kinase |
| Symbol: | Tek |
| UniProt ID: | Q02858 |
| Synonyms: | AA517024; Cd202b; Hyk; STK1; Tie-2; Tie2 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)