

Product datasheet for **GA204292**

Mouse **Tep1** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Tep1 CRISPRa kit - CRISPR gene activation of mouse telomerase associated protein 1
Symbol:	Tep1
Components:	GA204292G1, Tep1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA204292G2, Tep1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA204292G3, Tep1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P97499</u>
Synonyms:	Tp1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)