

Product datasheet for **GA204297**

Mouse Tes activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tes CRISPRa kit - CRISPR gene activation of mouse testis derived transcript
Symbol:	Tes
Components:	GA204297G1, Tes gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA204297G2, Tes gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA204297G3, Tes gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	D6Ert352e; Tes1; Tes2; TESS; testin2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)