

Product datasheet for GA204300

Mouse Prss40 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Prss40 CRISPRa kit – CRISPR gene activation of mouse protease, serine 40
Symbol:	Prss40
Components:	GA204300G1, Prss40 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074) GA204300G2, Prss40 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074) GA204300G3, Prss40 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074) 1 CRISPRa-Enhancer vector, SKU GEI00056 1 CRISPRa scramble vector, SKU GEI00077
UniProt ID:	<u>A6H6T1</u>
Synonyms:	Tesp2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

