

Product datasheet for **GA204303**

Mouse *Tex2* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tex2 CRISPRa kit - CRISPR gene activation of mouse testis expressed gene 2
Symbol:	Tex2
Components:	GA204303G1, Tex2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA204303G2, Tex2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA204303G3, Tex2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q6ZPJ0
Synonyms:	4930568E07Rik; AI553404; Def-5; Taz4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)