

Product datasheet for **GA204304**

Mouse **Tex261** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Tex261 CRISPRa kit - CRISPR gene activation of mouse testis expressed gene 261
Symbol:	Tex261
Components:	GA204304G1, Tex261 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA204304G2, Tex261 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA204304G3, Tex261 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q62302</u>
Synonyms:	3110001O07Rik; AA409339; AI480706; AL033351; TEG-261
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)