


Product datasheet for **GA204309**

Mouse Tex9 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tex9 CRISPRa kit - CRISPR gene activation of mouse testis expressed gene 9
Symbol:	Tex9
Components:	GA204309G1, Tex9 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA204309G2, Tex9 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA204309G3, Tex9 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9D845
Synonyms:	Gm19407; tsec-1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)