

Product datasheet for **GA204328**

Mouse Tgm3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tgm3 CRISPRa kit - CRISPR gene activation of mouse transglutaminase 3, E polypeptide
Symbol:	Tgm3
UniProt ID:	Q08189
Synonyms:	AI893889; TG(E); TGase-3; TGase E; TG E; TGE; we
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)