

## Product datasheet for GA204357

### Mouse Timp3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Timp3 CRISPRa kit - CRISPR gene activation of mouse tissue inhibitor of metalloproteinase 3
Symbol:	Timp3
UniProt ID:	<u>P39876</u>
Synonyms:	Timp-3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)