

Product datasheet for **GA204367**

Mouse Tle4 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | Tle4 CRISPRa kit - CRISPR gene activation of mouse transducin-like enhancer of split 4 |
| Symbol: | Tle4 |
| UniProt ID: | <u>Q62441</u> |
| Synonyms: | 5730411M05Rik; AA792082; Bce-1; Bce1; ESTM13; ESTM14; grg-4; Grg4; X83332; X83333 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)