

Product datasheet for **GA204373**

Mouse Otop1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Otop1 CRISPRa kit - CRISPR gene activation of mouse otopetrin 1
Symbol:	Otop1
Components:	GA204373G1, Otop1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA204373G2, Otop1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA204373G3, Otop1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q80VM9</u>
Synonyms:	A530025J20Rik; Otp1; tlt
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)