

## Product datasheet for **GA204377**

### Mouse Tspan7 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Tspan7 CRISPRa kit - CRISPR gene activation of mouse tetraspanin 7
Symbol:	Tspan7
UniProt ID:	<a href="#">Q62283</a>
Synonyms:	1200014P11Rik; A15; AI323365; Cd231; Mxs1; PE31; R74651; TALLA; Tm4sf2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)