

Product datasheet for **GA204379**

Mouse Tmod1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmod1 CRISPRa kit - CRISPR gene activation of mouse tropomodulin 1
Symbol:	Tmod1
UniProt ID:	P49813
Synonyms:	E-Tmod; Tmod
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)