

Product datasheet for **GA204383**

Mouse Tnnc1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tnnc1 CRISPRa kit - CRISPR gene activation of mouse troponin C, cardiac/slow skeletal
Symbol:	Tnnc1
UniProt ID:	P19123
Synonyms:	A1874626; cTnC; cTnl; TnC; tncc
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)