

## Product datasheet for **GA204394**

### Mouse Cd40 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Cd40 CRISPRa kit - CRISPR gene activation of mouse CD40 antigen
Symbol:	Cd40
UniProt ID:	<a href="#">P27512</a>
Synonyms:	AI326936; Bp50; GP39; HIGM1; IGM; IMD3; p50; T-BAM; Tnfrsf5; TRAP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)