

Product datasheet for **GA204401**

Mouse Cd40lg activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cd40lg CRISPRa kit - CRISPR gene activation of mouse CD40 ligand
Symbol:	Cd40lg
UniProt ID:	P27548
Synonyms:	CD40-L; Cd40l; CD154; gp39; HIGM1; IGM; IMD3; Ly-62; Ly62; T-BAM; Tnfsf5; Tnlg8b; TRAP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)